
MicroSet 40x-1024x

Art. No.XSP-42XT

- DE Bedienungsanleitung
- GB Operating Instructions
- FR Mode d'emploi
- NL Handleiding
- IT Istruzioni per l'uso
- ES Instrucciones de uso
- PT Manual de utilização

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RISK to your child.

Aids with sharp edges and tips are sometimes used with this device. Please store the device and all of its accessories and aids out of the reach of children. There is a risk of INJURY.

This device contains electronic components that are powered by either a mains connection or batteries. Never leave a child unsupervised with this device. The device should only be used as per these instructions otherwise there is a serious RISK of ELECTRICAL SHOCK.

Children should only use this device under supervision. Keep packaging materials (plastic bags, rubber bands, etc.) away from children. There is a risk of SUFFOCATION.

The chemicals and liquids provided should be kept out of reach of children. Do not drink the chemicals! Hands should be washed thoroughly under running water after use. In case of accidental contact with the eyes or mouth rinse with water. Seek medical treatment for ailments arising from contact with the chemical substances and take the chemicals with you to the doctor.

FIRE-/ DANGER OF EXPLOSION!

Do not expose the device to high temperatures. Use only the mains adapter supplied or those battery types recommended. Never short circuit the device or batteries or throw into a fire. Exposure to high temperatures or misuse of the device can lead to short circuits, fire or even explosion!

RISK of material damage

Never take the device apart. Please consult your dealer if there are any defects. The dealer will contact our service centre and send the device in for repair if needed.

Do not subject the device to temperatures exceeding 60° C.

TIPS on cleaning

Remove the device from its energy source before cleaning (remove plug from socket / remove batteries).

Clean the exterior of device with a dry cloth. Do not use cleaning fluids so as to avoid causing damage to electronic components.

Clean the lens (objective and eyepiece) only with the cloth supplied or some other soft lint-

free cloth (e.g. micro-fibre). Do not use excessive pressure - this may scratch the lens.

Dampen the cleaning cloth with a spectacle cleaning fluid and use it on very dirty lenses.

Protect the device from dust and moisture. Store the device in the bag supplied or in its original packaging. Batteries should be removed from the device if it is not going to be used for a long period of time.

DISPOSAL

Dispose of the packaging material/s as legally required. Consult the local authority on the matter if necessary.

Do not dispose of electrical equipment in your ordinary refuse. The European guideline 2002/96/EU on Electronic and Electrical Equipment Waste and relevant laws applying to it require such used equipment to be separately collected and recycled in an environment-friendly manner. Empty batteries and accumulators must be disposed of separately. Information on disposing of all such equipment made after 01 June 2006 can be obtained from your local authority.

Here are the parts of your microscope (Fig. 1-6):

- 1 10x WF Eyepiece
- 2 16x WF Eyepiece
- 3 Barlow Lens
- 4 MicroOcular
- 5 Eyepiece supports
- 6 Microscope Head
- 7 Set Screw
- 8 Objective Nosepiece
- 9 Objective
- 10 Clips
- 11 Microscope Stage
- 12 LED Illumination (transmitted light)
- 13 Microscope Base
- 14 Selection Knob for Illumination
- 15 Power supply
- 16 Focus knob
- 17 Color Filter
- 18 LED Illumination (reflected light)
- 19 Ulead PhotoExplorer SE Software
- 20 5 Slides, 5 Cover Sips and 5 Prepared Specimens plastic box
- 21 Specimens:
 - a) Yeast
 - b) "Gum Media" Glue
 - c) Sea Salt
 - d) Shrimp Eggs
- 22 Transformer with Plug and Power Cable
- 23 Specimen Slicer
- 24 Shrimp Hatchery

- 25 Tweezers
- 26 Pipette
- 27 Carrying Case

How do I use my microscope?

Before you assemble your microscope, make sure that the table, desk or whatever surface that you want to place it on is stable, and does not wobble.

Apart from that, you will need a power source for your microscope's light (220V-330V).

How do I operate the electric LED illumination?

For the illuminator, you need the transformer with the plug and power cable (Fig. 6, 22). First, place the small plug on the power cable into the power connection on the back of the microscope (Fig. 1, 15). Then, insert the large plug into a power outlet.

There are two lights on the microscope. They do not use light bulbs, but rather light-emitting diodes (LED). The first lamp shines onto the specimen from below and the second from above. (The thing that you want to observe with the microscope is called the object or specimen, by the way.) You can use each lamp on its own, or both of them together. There is a selection knob for this (Fig. 1, 14). It has three numbers: I, II and III.

If you select the ...

- I, the light only comes from below (transmitted light).
- II, the light only comes from above (reflected light).
- III, both lamps shine light on the specimen.

For transparent objects (transmitted-light objects), number I is best. In order to observe firm, non-transparent objects (direct-light objects), select number II. For semi-transparent objects, it is best to select number III.

It is not recommended to use number III for transmitted-light objects on slides, since the light may cause reflections on the surface of the slide, which will disturb your observation.

When do I use the color filters?

The color filters (Fig. 1, 17) are located below the microscope stage (Fig. 1, 11). They help you when you are observing very bright or clear specimens. Here, you can choose from various colors. This helps you better recognize the components of colorless or transparent objects (e.g. grains of starch, protozoa).

How do I adjust my microscope correctly?

First, loosen the screw (Fig. 1, 7) and turn the microscope head (Fig. 1, 6) into a comfortable viewing position.

Each observation starts with the lowest magnification.

Adjust the microscope stage (Fig. 1, 11) so that it goes all the way down to the lowest position. Then, turn the objective nosepiece (Fig. 1, 8) until it clicks into place at the lowest magnification (objective 4x).

Note:

Before you change the objective setting, always move the microscope stage (Fig. 1, 11) to its lowest position. This way, you can avoid causing any damage!

Now, insert the 10x eyepiece (Fig. 1, 1) into the Barlow lens (Fig. 1, 3). Make sure that the Barlow lens is placed all the way into the eyepiece supports (Fig. 1, 5) and is not pulled out (Fig. 2).

How do I observe the specimen?

After you have assembled the microscope with the adequate illumination and adjusted it

correctly, the following basic rules are to be observed:

Start with a simple observation at the lowest magnification. This way, it is easier to position the object in the middle (centering) and make the image sharp (focusing).

The higher the magnification, the more light you will require for a good image quality.

over the illumination (Fig. 1, 12).

In the next step, take a look through the eyepiece (Fig. 1, 1) and carefully turn the focus knob (Fig. 1, 16) until the image appears clear and sharp.

Now you can select a higher magnification by slowly removing the Barlow lens (Fig. 1, 3) from the eyepiece support (Fig. 1, 5). When the Barlow lens is almost completely pulled out, the magnification can be increased to almost double.

If you would like an even higher level of magnification, insert the 16x eyepiece (Fig. 1, 2)

From there, the light is reflected back and passes through the objective and eyepiece (where it gets magnified) into the eye. This is reflected light microscopy.

For transparent objections (e.g. protozoa), on the other hand, the light shines from below, through the opening in the microscope stage and then through the object.

The light travels further through the objective and eyepiece, where it is also magnified, and finally goes into the eye. This is transmitted-light microscopy.

Many microorganisms in water, many plant components and the smallest animal parts are already transparent in nature. Others have to be prepared. We may make them transparent through a treatment or penetration with the right materials (media), or by taking the thinnest slices from them (using our hand or a specimen slicer), and then examine them. You can read more about this in the following sections.

How do I make thin specimen slices?

Only do this with the supervision of your parents or another adult.

As I already pointed out, the thinnest slices possible are taken from an object. In order to get the best results, we need some wax or paraffin. It is best if you get a candle. Place

the wax in a pot and heat it carefully over a low burner. Now, dip the object in the liquid wax a few times. Then, let the wax get hard. Using the specimen slicer (Fig. 6, 23) or a knife/scalpel, cut the smallest slices from the object that is covered with wax. These slices are to be laid on a slide and covered with a cover slip.

How do I make my own specimens?

Take the object that you want to observe and place it on a glass slide (Fig. 6, 20). Then, add a few drops of distilled water on the object (Fig. 7) using a pipette (Fig. 7, 26). Now, place a cover slip vertically at the edge of the drop of water, so that the water runs along the edge of the cover slip. Then, slowly lower the cover slip over the water drops (Fig. 8).

Note:

The included glue "gum media" (Fig. 5, 21), is used to make permanent prepared specimens. Use this in place of the distilled water. If you want to keep the object in place on the slide permanently, use the gum media.

and turn the objective nosepiece (Fig. 1, 8) to a higher setting (10x or 40x).

Important tip:

The highest magnification is not always the best for every specimen!

Note:

Each time the magnification changes (eyepiece or objective change, pulling out the Barlow lens), the image sharpness must be readjusted with the focus knob (Fig. 1, 16). When doing this, make sure to be careful. If you move the microscope stage too quickly, the objective and the slide could come into contact and become damaged!

Which light for which specimen?

With this unit, a reflected light and transmitted light microscope, you can observe transparent, semi-transparent as well as non-transparent objects.

The image of the given object of observation is "transported" through the light. As a result, only the correct light will allow you to see something!

If you are observing non-transparent (opaque) objects (e.g. small animals, plant components, stones, coins, etc.) with this microscope, the light falls on the object that is being observed.

How do I install the MicroOcular?

Note:

The MicroOcular only works with the included Barlow lens! The magnification setting is changed when the MicroOcular is used, and it must be readjusted with the focus knob.

First, remove the Barlow lens (Fig. 4, 3) with the current eyepiece from the eyepiece support (Fig. 4) and place the MicroOcular (Fig. 5, 4) in the eyepiece (Fig. 5).

How do I install the software?

Do not plug the MicroOcular into your PC yet.

First, insert the included CD-ROM in your CD or DVD drive. The CD menu will open automatically. If the menu does not open, please open the Windows Explorer and select "CD Drive." Start the file "AutoRun.exe" by double clicking it.

Note:

If you already have an image-editing software program installed, you can skip to point 9.2. Only the driver for the MicroOcular will be installed. If you are not sure, continue with point 9.1.

How do I install the image-editing software Ulead PhotoExplorer 7.0 SE?

Please take the following steps:

1. Click on "Install PhotoExplorer 7.0 SE."
2. A choice of languages appears. Select the language that you prefer, and click „OK“ to confirm.
3. At "Welcome," click "Continue>". Now the License Agreement appears, which you must confirm with "Yes."
4. Now, please enter your name along with the serial number. You can leave the field "Company Name" blank if you wish. You'll find the serial number on the CD-ROM case. Confirm the information with "Continue>". If you have made a mistake while entering the serial number, a warning will appear. Click "OK" and enter the serial number correctly.
5. In the next 4 windows, you will be asked to select the „Target Path“, „Create Icons“, „Select Program Folder“, and „TV System.“ Confirm the information with "Continue>".
6. Now the "Setup Status" window appears.

The program files are being copied to your PC. This can take several minutes.

7. Once the window has disappeared, it may take a moment until the „File Allocation“ appears. Click "Continue>" if you would like to use the Ulead Software to view or edit all images on your PC in the future. If not, select „Do not Change Current File Allocation,“ and then click „Continue.“
8. If the Internet Explorer window now opens, please close it.
9. The window "InstallShield Wizard Complete" appears. Click on „No, I will restart the computer later“ and finally „Finish.“

The image editing software is now installed.

How does the driver installation for the MicroOcular work?

Make sure that the MicroOcular is not connected to your PC.

Please take the following steps:

1. Click on "Install Bresser PC-Ocular Driver."
2. A choice of languages appears. Select the language that you prefer, and click „OK“ to confirm.
3. At "Welcome," click "Continue>".
4. Select "Complete" and confirm with "Continue>".

5. In the "Start Copying," select "Continue>" again.
6. Now the "Setup Status" window appears. The program files are being copied to your PC. This can take several minutes.
7. Once the window has disappeared, it may take a moment until the „InstallShield Wizard Complete“ appears. Click on "Finish."

The driver is now installed. Click on "Exit" in the menu.

Note:

In order to work with the MicroOcular for the long-term, I recommend that you always connect the MicroOcular to the same USB port on your PC.

How do I assign the driver?

Restart your PC. After rebooting the system, connect the MicroOcular to a free USB port.

The MicroOcular will now be assigned the driver. In Windows XP, a balloon appears at the lower right of the screen. The display will switch between "CIF Single Chip" to "MicroOcular" for a short period of time. Then the balloon will read "The hardware is now installed and is ready to use."

The MicroOcular is now ready to use.

How do I use the MicroOcular?
Preparation

1. First, place a specimen in the microscope and focus the image.
2. Remove the eyepiece and the Barlow lens from the eyepiece support, and remove the protective cover from the MicroOcular and place it in the eyepiece support.
3. Start your PC, if it is not already on, and connect the MicroOcular to your PC's USB port.

Presenting and saving images from the MicroOcular on your PC

1. Start Ulead PhotoExplorer 7.0 SE.
2. Select "Video" from the menu list and then "Record Picture."
3. The window "Record Picture from Video Unit" appears (Fig. 9). Click "Live View." Now, the image of the MicroOcular should be visible on the screen. If not, click "Setup" and then "Options." The "Recording Unit" window appears. In "Select Camera," choose „Microocular“ and confirm with „OK.“
4. Focus the image once more. As soon as the desired image is visible, click "Record Now."
5. The image is now saved to your computer. Click on "Close" and the picture appears in Ulead PhotoExplorer.

Fig. 9

Ulead PhotoExplorer 7.0

You'll find a detailed instruction manual in PDF format for Ulead PhotoExplorer 7.0 SE on the CD in the folder "Document\Manual." In this folder, there are folders for a variety of languages. Should your language not be available, make a direct request to Ulead to see if they have a version in that language available.

Troubleshooting

Error	Solution
No recognizable image	• Turn on light • Readjust focus
Image flickers (while observing with MicrOcular) = Monitor refresh rate not adequate	• If necessary, increase resolution of graphics card
Software installation reports "no XP approved"	• Confirm with "OK"

Make sure your microscope has a long service life.

Clean the lens (objective and eyepiece) only with the cloth supplied or some other soft lint-free cloth (e.g. microfibre). Do not press hard as this might scratch the lens.

Ask your parents to help if your microscope is really very dirty. The cleaning cloth should be moistened with cleaning fluid and the lens wiped clean using little pressure.

Make sure your microscope is always protected against dust and dirt. After use leave it in a warm room to dry off. Then install the dust caps and keep it in the case provided.

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