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


ENGLISH

We appreciate your choice on our Konus Microscope.

This leaflet will explain how to install this instrument and what technical characteristics it has.

TECHNICAL CHARACTERISTICS

- Stereoscopic microscope with eyepieces inclined at 45°.
- Pair of W10x eyepieces
- Two couples of achromatic objectives (2x and 4x) on rotating turret.
- Incidental or transmitted light with bulbs 12V/10W.
- External transformer for energy supply.
- Two specimen-carrying discs (black/white and transparent)
- Rack focus.
- Rubber eyecups on the eyepieces.
- Dioptic regulation device.

COMPONENTS (see figure):

- 1) Couple of eyepieces
- 2) Head
- 3) Stand
- 4) Base
- 5) turret objective holder.
- 6) 2 couples of objectives.
- 7) Specimen disk (1 b/w and 1 transparent)
- 8) (internal) bulb for transmitted light.
- 9) Bulb for incidental light.
- 10) Eyepiece holder.
- 11) 3 Position switch.
- 12) Rack and pinion focus.
- 13) Rubber eyecups on the eyepieces.
- 14) Dioptic regulation device.

MOUNTING

- Open the packaging and take out all the components.

Position the base (4) with the stand (3) and head (2) on a dry and stable table.

-Take the transformer and plug it into an electric socket, (make sure that the data on the label of the transformer matches that of your electrical supply).

-Insert the jack plug of the transformer into the socket of the base (4) of the instrument.

- Turn on the light with the three way switch (11); with this switch, it is possible to turn on the two bulbs separately or together. Now the microscope is ready for use.

NOTE.

Before proceeding with your observation, it is necessary to regulate the microscope with the following operations.

INTER-PUPILLARY DISTANCE REGULATION: For a better field of vision, it is impor-

STEREOSCOPISCHE MICROSCOOP KONUS DIAMOND

tant that the distance between the two eyepieces (1) are the same as your eyes. To do this, (the microscope must be on, with the eyepieces mounted) you must look into the eyepieces, if you cannot see a single image, you must adjust the two holders of the eyepieces.

DIOPTRIC REGULATION.

One of the two eyepiece holders (14), has the possibility of turning around, to higher or lower the eyepiece. This function permits a perfect vision even for those that present dioptic differences from one eye to the other. For a perfect regulation, take the following steps:

1. Looking into the eyepiece that cannot be trimmed, bring in the image of the specimen using only the focusing knob, macrometric (12).

2. Look into the trimmable eyepiece (14), with the other eye, and bring into focus the specimen turning the eyepiece holder.

N.B: A good preliminary regulation of the instrument (Dioptic and distance regulation) is necessary for a clear and easy vision. If more people use this instrument, it will be necessary to use this process each time.

ENLARGEMENTS: The Konus Diamond model is characterised of having a rotatable objective turret with two pairs of achromatic objectives. To select one pair of objectives or the other, the turret need only be turned until a click is heard.

To calculate the enlargements use the formula:

$$\text{ENLARGEMENTS} = \frac{\text{MAGNIF. OF EYEPIECES} \times \text{MAGNIF. OF OBJECTIVE}}{1}$$

In the Konus Diamond model with 10x eyepieces and 4x objectives 10x2 and 10x4 will give us enlargements of 20 times and 40 times the original. It is possible to vary the enlargements changing the ocular lens or the objectives (the objectives are those with a dovetail finish and are easily interchangeable).

MAINTENANCE: The instrument does not present itself with particular maintenance problems, it is important to follow these simple rules:

Do not tamper with the instrument, especially the electrical wiring and the optical parts.

Protect the instrument from dust, acids, shocks and/or impacts, water and/or humidity.

Do not place the transformer within reach of children.

The eventual vision by children of the microscope should be taken under the supervision of adults.

OPTIONAL ACCESSORIES:

- Pair of eyepieces 5x #5470
- Pair of eyepieces 15x #5463
- Pair of eyepieces 20x #5433
- Photography adaptor #5124 (T2 ring not included)
- CCD camera adaptor #1599
- CCD camera #5821
- Dark field set #5138.

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© KONUS 2012 ALL RIGHT RESERVED - PRINTED IN P.R.C.
KONUS IS REGISTERED TRADE MARK OF KONUS IN U.S.A., ITALY, SPAIN AND